The natural way of Sweden

SISSEL® SPINEFIT

- Gebrauchsanweisung
- Instructions for use
- Mode d'emploi
- Gebruiksaanwijzing
- Instrucciones de uso
- Istruzioni per l'uso
- Bruksanvisning

4 250694 714127

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
The SISSEL® SpineFit is an innovative fitness tool combining both muscle release and effective training

- intensive, deep penetrating pressure on the entire musculature alongside the spinal column
- supports mobility of the joints and muscle relaxation
- noticeable relaxation of the shoulder girdle already in resting position
- developed based on everyday practical experience in therapy & studio by Katrin Zinke (Certified Physical Education Teacher & Shiatsu Therapist) and Nina Metternich (Physiotherapist & Osteopath)
- unique sports tool that can be used in physiotherapy, personal or group training

The number, arrangement and weight of the balls created an effective training device that stimulates the energy points alongside the spinal column, bringing the body back into balance.

Katrin Zinke

Nina Metternich

SISSEL® SpineFit – Noticeable back mobilization with an intensive feel-good effect

In our daily work in therapy and group fitness, we use balls of different sizes and firmness to support exercises. With the SISSEL® SpineFit we found a solution to expand and ideally support therapy and training by using balls in an optimal arrangement, size, firmness and weight.

The unique shape of the SISSEL® SpineFit with its ease of use improves the shoulder-neck area. Mechanical stimuli, initiated by movements on the SISSEL® SpineFit, open the small vertebral joints of the spinal column and, thereby, can release blockages.

With its ease of use, the unique shape of the SISSEL® SpineFit provides relaxation in the shoulder-neck area. Resting in supine position on the SISSEL® SpineFit has a stimulating effect on the energy points of the body, thus supporting the energetic balance. The SISSEL® SpineFit is a versatile tool for therapy, personal and group training.


It is a useful addition to back, Pilates and yoga exercises. By facilitating joint alignment and challenging positions (balance), new stimuli are set. The intuitive use and quickly noticeable effect make the SISSEL® SpineFit an optimal sports tool for regular home training.

Exercise video at youtube.com/sisselcom

Exercise options

Relaxing position

In this position, you relax the muscles alongside the spine.

Execution

- lying in neutral supine position, lengthwise on the SISSEL® SpineFit

Bridging

Mobilizes the small vertebral joints and makes the acupressure effect on the energy points noticeable.

Execution

- spine rests lengthwise on the SISSEL® SpineFit, feet hip-width on the mat
- roll up vertebra by vertebra up to shoulder blades and down again, perceiving the mobility of the spine

Shoulder circles

Relieves tensions in the shoulder and neck muscles.

Execution

- spine rests lengthwise on the SISSEL® SpineFit, hands pointing to the ceiling, elbows slightly bent
- make circular movements out of the shoulder girdle with straight arms; thereby, you perceive the mobility of the shoulder blade with the release of tense muscles

Exercise options

Arm circles

Mobilizes the shoulder joints and relieves tensions in the capsule.

Execution

- spine rests lengthwise on the SISSEL® SpineFit, hands pointing to the ceiling, elbows slightly bent
- circle your arms in a wide arc and change direction

Massage of the shoulder girdle

Stimulates trigger points (energy points) and relieves tensions in the shoulder girdle.

Execution

- supine position with the SISSEL® SpineFit crosswise underneath the top of the scapula, legs set up, upper body slightly raised, head in your hands
- slowly roll from the outside of the scapula to the opposite side and back

Package


Stimulates energy points of the sacrum and releases tension of the deep vascular musculature.

Execution

- supine position with the SISSEL® SpineFit diagonally underneath the sacrum, lower legs hang above the pelvis, hands on the SISSEL® SpineFit stabilizing the upper body
- circle knees to mobilize the gluteal muscles and the sacrum, change direction

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SISSEL® GmbH